

THE 522KL

DIGITAL PEDIATRIC TRAY SCALE

The 522KL pediatric scale was specifically created for nurseries, delivery rooms, and pediatric practices, making it an ideal choice for newborns, infants, and toddlers. With a maximum capacity of 50 pounds, this scale offers providers and patients a range of benefits, including safety, comfort, and flexibility. The scale ensures a secure platform, keeps patients calm for accurate readings.

The 522KL display is positioned 17" above the scale base, ensuring an unobstructed easy-to-read weight measurement. The user-friendly digital interface includes features such as Zero, Tare, Pre-Tare, Hold/Release, and Auto-off, making it convenient and easy to use. Its motion detection technology ensures accurate weight measurements even when the patient moves.

The 522KL provides easy integration and seamless communication with today's EMR applications with the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices.

The 522KL pediatric scale offers a built-in measuring tape, allowing users to measure weight and height simultaneously. It also provides the flexibility of choosing between an optional digital or mechanical height rod.

522KL-EHR

Scale with attachable digital height rod



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 522KL is your choice for safety, comfort, and reliability.

120+ LOCAL
FIELD SPECIALISTS

30-60 SEC
AUTO WEIGHT HOLD

#1 TRUSTED
BRAND

100 YEARS
OF SERVICE

Health o meter®
Professional

Health o meter® Professional 522KL

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

50 lb / 23 kg

Resolution

0 - 20 lb / 0.2 oz; 20 - 50 lb / 0.5 oz;
0 - 9 kg / 5g; 9 - 23 kg / 10 g

EMR Connectivity

Via optional Pelstar®
wireless technology

Display

1" / 25 mm, LCD display

Power Source

6 AA batteries included,
100-240V adapter (ADPT31)
not included

Functions

LB / KG conversion
LB / KG lock out
Zero
Tare
Pre-Tare
Hold / release
Reweigh
Auto zero
Auto off
Motion detection
Audible/mute sound option
Variable auto-off time
Enabling unit view with unit of measure locked
Auto-weight hold for 30 to 60 seconds
Everlock® - permanent weight lock

Tray Size (w x d x h)

24 1/8" x 14 1/2" x 2 5/8"
(613 mm x 368 mm x 67 mm)

Product Footprint

24 1/8" x 21" x 23 1/8"
(613 mm x 533 mm x 587 mm)

Product Weight

20 lb / 9 kg

Measuring Tape

0" - 23" / 0 cm - 58 cm

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
522KL	Digital Pediatric Tray Scale, LB & KG
522KG	Digital Pediatric Tray Scale, KG only
522KL-HR	522KL with attachable mechanical height rod
522KG-HR	522KLG with attachable mechanical height rod
522KL-EHR	522KL with attachable digital height rod
522KG-EHR	522KG with attachable digital height rod
522HR	Mechanical height rod for 522 series of scales
522EHR	Digital height rod for 522 series of scales
SS-522KL	ScaleSurance® 2 year extended warranty for 522 series*
ADPT31	Optional adapter for 522 series of scales

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscales.com/ConnectivitySolutions



To learn more, call **1-800-815-6615**
or visit homscales.com/522KL

Health o meter®
Professional



*Not available in all countries.