

750 Garcia Ave, Pittsburg CA 94565 +1.925.427.7800 HSIheadwalls.com

Cogent[™] Slimline

Training Headwall

The Non-Functional Slimline Training Headwall is designed to give the aesthetic of a patient room without the cost associated with medical gas simulation.





Features/Specifications:

- These units can be wired for Electrical capabilities [inquire at time of order]
- The Non-Functional Slimline Training Headwall can be moved with ease when lab spacing needs change
- The Non-Functional Slimline Training Headwall horizontal rail is 52" wide x 4" Deep x 10" High
- Installation is simple: locate the Non-Functional Slimline Training Headwall in position and secure it to the building wall.

Options:

Add:

- Overbed Light Fixture [with pull cord switch]
- Integrated Non-Functional Nurse Call*
- Integrated Non-Functional Code Blue*
- · Vacuum Regulator
- · Vacuum canister and holder