

# ENT Chair with Foot Operated Pump

## Model 8612


The UMF Medical 8612 offers a manually adjustable backrest, footrest, and foot operated pump for hi-lo movement. Offers a 360 degrees rotation, hinged arms, and 250lbs weight capacity.



### Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



River Rock



Sand Gray



Steel Blue



Midnight Blue

### Key features at a glance

- › 250lb (113.4kg) patient weight capacity
- › Ultra comfort seamless top with antimicrobial guard
- › Sturdy steel frame
- › Available in all UMF Medical upholstery options
- › 360° rotation with swivel lock
- › 10° elevation from horizontal position
- › Manual foot pedal height adjustment
- › Manually adjustable backrest
- › Manually adjustable footrest
- › Hinged arms
- › 3-year limited warranty

### Specifications

Min. Seat Height: 23" (58.4cm)  
 Max. Seat Height: 38" (96.5cm)  
 Overall Width: 31.5" (80cm)  
 Seat Width: 18" (45.72cm)  
 Arm Length: 23" (58.4cm)  
 Length: 71" (180.3cm) fully extended  
 Net Weight: 166lb (75.3kg)  
 Shipping Weight: 225lb (102kg)  
 Shipping Dimensions: L44xW34xH64  
 (111.8cm x 86.4cm x 162.6cm)

### Accessories

031 – Optional backrest with articulating headrest