

Phlebotomy Chair with Foot Operated Pump

Model 8675


The UMF Medical's 8675 uniquely offers an interlinked backrest and leg section for reclining patient to a horizontal position; as well as, manual foot pedal height adjustment. The 8675 has a 360 degree rotation with swivel lock, and 250lbs weight capacity.



Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



River Rock



Sand Gray



Steel Blue



Midnight Blue

Key features at a glance

- › 250lb (113.4kg) patient weight capacity
- › Ultra comfort seamless top with antimicrobial guard
- › Available in all UMF Medical upholstery options
- › Sturdy steel frame
- › 360° rotation with swivel lock
- › 10° elevation from horizontal position
- › Manual foot pedal height adjustment
- › Interlinked backrest and leg section for reclining patient to a horizontal position
- › Hinged arms
- › 3-year limited warranty

Specifications

Length: 67" (170.2cm) fully extended
 Min. Seat Height: 23" (58.4cm)
 Max. Seat Height: 38" (96.5cm)
 Overall Width: 31.5" (80cm)
 Seat Width: 18" (45.72cm)

Arm Length: 23" (58.4cm)
 Net Weight: 171lb (77.6kg)
 Shipping Weight: 230lb (104.4kg)
 Shipping Dimensions: L44xW34xH64
 (111.8cm x 86.4cm x 162.6cm)