

## Family Practice Table

### MODEL 8870

- Premium, vinyl upholstery over 2" firm, foam padding (5 cm)
- Full 28.5" width for patient comfort (72.39 cm)
- Air-spring adjustable backrest
- Coated, two position adjustable stirrups
- Pull-out leg rest on steel track
- Easy care, Gray, laminate base
- Color coordinated doors and top available
- Chip resistant edges
- 1 large side storage compartment with shelf
- Storage compartment is accessible from both sides (no right or left handed tables)
- 2 drawers with all steel slides
- Built with plywood and laminate panels that are California CARB II compliant
- Top sections are easily replaceable in the field
- Adjustable feet for leveling
- Paper dispenser included
- 400 lbs. load capacity under normal use (181.43 kg)



#### SPECIFICATIONS

58"/72" extended L x 28.5" W x 31" H



California residents please be advised of the following, pursuant to Proposition 65:

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

### Features




Pull-out, slate, laminate, leg rest



2 drawer storage



2-Position Stirrups



Paper Roll Storage



Airlift Backrest




Both Sides Access

### Laminate Color



Gray

### Option



040- Paper Holder

### Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black