


DISPOSABLE PENLIGHTS

ADLITE™

- Sealed concentric beam, disposable diagnostic penlight
- Activated by depressing pocket clip
- Available with or without pupil gauge
- Measures 4 1/2" long, 1/2" diameter
- Sold in 6 pack



ADLITE™

Diagnostic Penlights

ORDERING INFORMATION

#351	ADLITE™ Disposable Penlight, 6/Pk
#351P	ADLITE™ Disposable Penlight, w/Pupil Gauge, 6/Pk

ADLITE™ PLUS

- Sealed concentric beam, disposable diagnostic penlight
- Activated by depressing pocket clip
- White available with or without pupil gauge
- Measures 4 1/2" long, 1/2" diameter
- Sold individually or in display cylinder of 22
- Available in 8 colors




ORDERING INFORMATION

#356*	ADLITE™ Disp. Penlight, Specify Color
#358AST	ADLITE™ Disp. Penlight, Assorted
#358TVB	ADLITE™ Disp. Penlight, BK, Purp, Teal
#358W	ADLITE™ Disp. Penlight, White
#358WP	ADLITE™ Disp. Penlight, White/Pupil Gauge

(358AST, 358TVB, 358W, 358WP, All come 22 pieces in cylinder sleeve)

Color Choices



ADC
AMERICAN DIAGNOSTIC CORPORATION
55 Commerce Drive
Hauppauge, New York 11788
tel: 1-800-232-2670
fax: 631-273-9659

www.adctoday.com



REUSABLE PENLIGHTS

Four distinctly different models to suit your personal preferences.



ADLITE PRO™

Our finest pocket light, and physician's choice for evaluation of the pupillary reflex.

- LED lamp for bright, white illumination for true tissue rendering
- LED lamp provides 10,000 hour use
- Durable CNC machined brass housing in a striking black matte finish
- Positive feel power switch
- Measures 5 5/8" long, 1/2" diameter
- Convenient Pocket Clip
- Complete with 2 "AAA" batteries and an attractive gift box
- Lifetime warranty on lamp cartridge, 2 year warranty on penlight

ORDERING INFORMATION

#355BK	Adlite Pro™ Physicians Penlight
--------	---------------------------------



METALITE II™

- Brushed aluminum housing
- Plunger activated
- Measures 5 1/4" long by 1/2" diameter
- Two "AAA" batteries included

ORDERING INFORMATION

#353BK	Metalite II™ Plunger Penlight, Black
#353RB	Metalite II™ Plunger Penlight, Royal Blue
#353R	Metalite II™ Plunger Penlight, Red
#353S	Metalite II™ Plunger Penlight, Silver



ADLITE II™

- Plastic housing for lightweight comfort
- Pocket clip activated
- Measures 5" long by 1/2" diameter
- Two "AAA" batteries included

ORDERING INFORMATION

#354	ADLITE II™ Reusable Penlight
#354CYL	ADLITE II™ Reusable Penlight, 18/cylinder



METALITE™

- Brushed aluminum housing
- Locking push toggle activated
- Measures 5 3/4" long by 1/2" diameter
- Two "AAA" batteries included



ORDERING INFORMATION

#352BK	METALITE™ Penlight, Black
#352W	METALITE™ Penlight, White
#352WP	METALITE™ Penlight, White/Pupil Gauge
#352	METALITE™ Penlight, Silver
#352CYL	METALITE™ Penlight, 22/Cylinder (Silver)
#352AST	METALITE™ Penlight, 22/Cylinder (Assort colors)

Display packaging available. Add "Q" suffix to order number.