



### P-7202

Cardiac & Anesthesia Cabinet

**Overall Size**

28"W x 18"D x 33"H

**Top Drawer**

21¼"W x 16"D x 3¼"H

**Bottom Drawer**

21¼"W x 16"D x 3¼"H



*shown with optional cart accessories*

### P-7202-C

Same as P-7202 w/ complete accessory package

**optional accessories for P-7202**

- O<sub>2</sub> holder (holds 2 E-size tanks)
- 25' electrical outlet cord
- I.V. rod
- Push handle
- Cabinet locking rod



## anesthesia & treatment cabinets

### P-7202 features

- Designed for cardiac and other emergency uses or as a large anesthesia cabinet
- Provides space for instruments, vials, and other necessary equipment
- 4" double ball bearing swivel casters



**new!**  
co-re  
caster™

### P-1100-SS

Anesthesia Table with Drawer

**Overall Size**

21<sup>3</sup>/<sub>8</sub>"W x 14"D x 34<sup>1</sup>/<sub>4</sub>"H

**Trays**

12<sup>5</sup>/<sub>8</sub>" W x 5<sup>5</sup>/<sub>8</sub>"D x 19<sup>5</sup>/<sub>8</sub>"H

**Drawer**

15<sup>3</sup>/<sub>8</sub>" W x 12<sup>1</sup>/<sub>2</sub>" D x 3<sup>1</sup>/<sub>4</sub>"H

**features**

- Ideal cart for storage of various instruments and supplies.
- Two removable stainless steel trays
- Towel bar and two mask hooks
- 2" double ball bearing swivel casters



### P-5090

Small Treatment Cabinet

**Overall Size**

16"W x 20"D x 33"H

**Drawer**

13<sup>1</sup>/<sub>2</sub>"W x 12"D x 3<sup>1</sup>/<sub>2</sub>"H

**features**

- Solid stainless steel retaining rail
- Internal shelf for optimal storage
- Optional locking rod
- 2" double ball bearing swivel casters



**new!**  
co-re  
caster™

### P-1105-SS

Anesthesia Cabinet with Drawers

**Overall Size**

19<sup>1</sup>/<sub>8</sub>"W x 15<sup>5</sup>/<sub>8</sub>"D x 33<sup>3</sup>/<sub>8</sub>"H

**Upper 3 Drawers**

13<sup>1</sup>/<sub>2</sub>"W x 12"D x 3<sup>1</sup>/<sub>2</sub>"H

**Bottom Drawer**

13<sup>1</sup>/<sub>2</sub>"W x 12"D x 7"H

**features**

- Provides extra storage space
- Full extending drawer hangers
- Removable 12<sup>5</sup>/<sub>8</sub>" x 19<sup>1</sup>/<sub>8</sub>" stainless steel tray
- 3" double ball bearing swivel casters



### P-5091

Large Treatment Cabinet

**Overall Size**

28"W x 18"D x 33"H

**Top Drawers**

21<sup>1</sup>/<sub>4</sub>"W x 16"D x 3<sup>1</sup>/<sub>4</sub>"H

**Bottom Drawer**

21<sup>1</sup>/<sub>4</sub>"W x 16"D x 5<sup>3</sup>/<sub>8</sub>"H

**features**

- Solid stainless steel retaining rail
- Internal shelf for optimal storage
- Optional locking rod
- 2" double ball bearing swivel casters