

sterile processing equipment



shown with hospital accessories

CDS-3260

Central Supply work table

Overall size
60"W x 32"D x 34"H

features

- Stainless steel construction
- Height adjustable floor glides



shown with hospital accessories

CDS-160

Instrument Container Wash Cart

Overall size
49"W x 26"D x 67"H

features

- Stainless steel construction
- Fully cart-washable
- Built to Steril-Gard specifications
- Fully welded frame
- Donut bumpers on legs
- 5" diameter stainless steel swivel casters
- Five stainless steel wire shelves
- Six adjustable container lid dividers



Caster with brake

Optional accessories

- 4" diameter swivel casters with brakes
- Sterile Wrap hangers measure 56"W x 30½"H
- Full-length lower shelf
- Tote box drawers:
5" or 10" deep
- Stainless steel drawers:
15½"W x 12½"D x 3½"H (installed, complete)