



CDS-153 Cart

Overall Size

36½"W x 26¾"D x 35½"H

- Optional linen rings to hold hamper or trash bags
- Linen ring locks in up position and swings down when not in use
- Equipped with rails on the all three surfaces
- Adjustable stainless steel center shelf
- 5"- diameter chrome plated swivel casters

multi-purpose carts

features

- All welded stainless steel construction for one single piece frame
- 5" chrome plated casters with diagonally mounted swivel locks provide stability
- Wrap around bumper mounted away from frame for easy cleaning
- Optional shelves and accessories available to configure the cart to suit individual needs



Swivel Plate Caster



Cart shown with optional ring