

P-5182 Dimensions

Standard Top 34"W x 42"D

Pedigo's Over-Operating Specialty Table is constructed of stainless steel and features easy, one-hand height adjustment.

- 3" total lock casters, with brakes and swivel locks
- Adjustable height from 37¹/₄" to 52¹/₄"high from floor
- 27¾" clearance between legs
- Other sizes available upon request.



P-5190-E Dimensions

Standard Top 36"W x 42"D

over-operating

la	\cup	l	U 5)

Dimension

28" ID Clearance

34" ID Clearance

P-5182-36	Base ONLY (X-Wide)	36" ID Clearance						
			-					
			Base (P-5182-3)	Base (P-5182-34)	Base (P-5182-36)			
P-5182-2	Top ONLY (Standard)	34"W x 42"D	•					
P-5182-3420	Top ONLY	34"W x 20"D	•					
P-5182-3424	Top ONLY	34"W x 24"D	•					
P-5182-3430	Top ONLY	34"W x 30"D	•					
P-5182-3436	Top ONLY	34"W x 36"D	•					
P-5182-3448	Top ONLY	34"W x 48"D	•					
			•					
P-5182-3820	Top ONLY	38"W x 20"D	•					
P-5182-3824	Top ONLY	38"W x 24"D	•					
P-5182-3830	Top ONLY	38"W x 30"D	•					
P-5182-3836	Top ONLY	38"W x 36"D	•					
P-5182-3842	Top ONLY	38"W x 42"D	•					
P-5182-3848	Top ONLY	38"W x 48"D	•					
			•					
P-5182-4220	Top ONLY	42"W x 20"D	•	•	•			
P-5182-4224	Top ONLY	42"W x 24"D	•	•	•			
P-5182-4230	Top ONLY	42"W x 30"D	•	•	•			
P-5182-4236	Top ONLY	42"W x 36"D	•	•	•			
P-5182-4242	Top ONLY	42"W x 42"D	•	•	•			
P-5182-4248	Top ONLY	42"W x 48"D	•	•	•			
P-5182-4620	Top ONLY	46"W x 20"D	•	•	•			
P-5182-4624	Top ONLY	46"W x 24"D	•	•	•			

The various sizes/configuration shown in the chart above are **ONLY available for the P-5182** over-operating table.

46"W x 30"D

46"W x 36"D

46"W x 42"D

46"W x 48"D

Adjustable Over-Operating Table may be height adjusted during the surgical case while the staff is gowned and sterile. To raise or lower the main table surface, simply engage the electric motor by depressing the foot pedal.

P-5182-4630

P-5182-4636

P-5182-4642

P-5182-4648

Top ONLY

Top ONLY

Top ONLY

Top ONLY

Part Number

P-5182-3

P-5182-34

Description

Base ONLY (Standard)

Base ONLY (Wide)

- 3" total lock casters, with brakes and swivel locks
- 28⁷/₈" clearance between legs
- Adjustable height from 40½" to 51½" high from floor (underside of table height from 35½" to 46½" high from floor)