

**P-1080-CV**

Height adjustable from 71" to 124"

**P-1080-6**

Height adjustable from 72½" to 124"



**P-1083-CV**

Height adjustable from 71" to 124"



**P-1083-6**

Height adjustable from 72½" to 124"

# infusion pump stand

**common features**

- Bag Drop Prevention Device allows for controlled descent of pole.
- Heavy duty 3" (8cm) sealed casters with precision ball bearings.
- Extremely durable powder coat finish base.
- Low center of gravity design and 27 lbs. base offers improved stability
- Pole sleeve welded securely to the base for reinforced stability.
- Small 23" (58cm) base reduces stumbles and falls while patient ambulates.
- Hand operated easygrip™ knob.

**P-1080-CV**

- High impact molded polymer IV hook assembly features six hooks.

**P-1080-6**

- Six stainless steel IV hooks

**P-1083-CV**

- High impact molded polymer IV hook assembly features six hooks.
- Pole made from copper nickel alloy.

**P-1083-6**

- Pole made from copper nickel alloy.
- Six stainless steel IV hooks



### Pole


Stainless steel resists staining, corrosion and rusting.

### Stainless Steel Base

Low center of gravity design and 27 lbs. base offers improved stability.

### 23" Diameter Base

Small base reduces stumbles and falls while patient ambulates.



Height adjustable from 70½" to 124"

## P-1080-SS

infusion pump stand

### features

- Bag Drop Prevention Device allows for controlled descent of pole.
- Stainless steel IV hook assembly features six hooks.
- Heavy duty 3" (8cm) sealed casters with precision ball bearings.
- Corrosion resistant all stainless steel base.
- Pole sleeve welded securely to the base for reinforced stability.



Stainless base

# infusion pump accessories



### PATIENT HANDLE

Adjustable handle provides easy mobility.


- P-3305 Stainless Steel
- P-3306 Copper nickel alloy
- P-3307 Black
- P-3305-T Tray Insert



### CV HOOK

Molded polymer six hook IV assembly

- P-1090-CV-TH Assembly



### INFUSION PUMP HOLDER

Holds up to 3 pumps

- P-1080-G Infusion Pump Holder



### 8-HOOK ASSEMBLY

Extended IV Hook Assembly

- P-1080-H Stainless Steel



### IV SET ORGANIZER

Also holds secondary bags

- P-1080-B Stainless Steel



### S'MORE HOOK

For mounting additional IV bags

- P-1080-J Stainless Steel

# infusion pump accessories



## WIRE BASKET

P-3201 Large (Stainless)  
P-3202 Large (Black)  
P-3211 Med (Stainless)  
P-3212 Med Black



## FOLEY HOOK

P-1080-D Foley 2 Hooks  
P-1080-E Foley 1 Hook  
*(Factory installed)*



## O<sub>2</sub> HOLDER

Holds E-size oxygen tanks  
P-3601 Stainless Steel



## OUTLET STRIP

Hospital grade  
P-3401 Six outlet strip



## INFUSION PUMP HOOP

Holds up to 6 pumps  
P-1080-F Stainless Steel  
P-1083-F Copper nickel alloy



## MULTIPLE PUMP HOLDER

Infusion pump holder  
P-1080-C Stainless Steel  
P-1083-C Copper nickel alloy

