



shown with optional
bassinet and pad

P-1110-B-SS
bassinet stand with drawer

new!
**co-re
caster™**

bassinet

P-1110-B-SS features

- Side mounted drawer measures:
16½"W x 21¼"L x 3⅛"H
- Design allows close approach from both sides
- Basket can be raised at either end
- 3" double ball bearing swivel casters

frame size

18¼"W x 31"L x 37¾"H
(Height with bassinet basket: 41"H)



shown with optional
bassinet and pad

P-1110-SS
bassinet stand

features

- Design allows close approach from both sides
- Basket can be raised at either end
- 3" double ball bearing swivel casters

frame size

- 18¼"W x 31"L x 37¾"H
- Height with bassinet basket: 41"H



shown with optional
bassinet and pad

P-1110-A-SS
bassinet stand with drawer

features

- End mounted drawer measures:
16½"W x 13¼"L x 7⅜"H
- Design allows close approach from both sides
- Basket can be raised at either end
- 3" double ball bearing swivel casters

frame size

- 18¼"W x 31"L x 37¾"H
- Height with bassinet basket: 41"H



shown with optional bassinet and pad

P-1115-A-SS

bassinet cabinet

features

- Bassinet carriage slides back to allow working area
- Friction lock secures the bassinet in both open and closed positions
- Removable internal shelf, positioned in middle left compartment
- 4" double ball bearing swivel casters

dimensions

19³/₄"W x 30⁷/₈"L x 37¹/₄"H

- Height with bassinet 44¹/₈"

cabinet dimensions

- 17"W x 29¹/₂"L x 24⁷/₈"H

drawer frame with removable pan

- 11¹/₄"W x 16¹/₂"L x 4"H



shown with optional bassinet and pad

P-1115-SS

bassinet cabinet

features

- Bassinet carriage slides back to allow working area
- Friction lock secures the bassinet in both open and closed positions
- Removable internal shelf, positioned in middle left compartment
- 4" double ball bearing swivel casters

dimensions

19³/₄"W x 30⁷/₈"L x 37¹/₄"H

- Height with bassinet 44¹/₈"

cabinet dimensions

- 17"W x 29¹/₂"L x 24⁷/₈"H

drawer frame with removable pan

- 11¹/₄"W x 16¹/₂"L x 4"H



ST-1833-SS

pediatric scale table

features

- Designed for use with common electronic scales
- Oversize 4" casters for easy movement
- One brake caster for stability

dimensions

33"W x 18"D x 35¹/₂"H

OPTIONS



Part No. 53
Bassinet Basket

basket measures

16"W x 30¹/₄"L x 9⁵/₈"H



Part No. 54
Bassinet Pad - Grey

pad measures

12¹/₂"W x 26¹/₂"L x 1" Thick