

**P-1015-TC**

**Dimensions**

28"D x 14¾"W x 40¼"H

- Stainless Steel
- Holds up to six P-1015
- 3" Genesis diameter swivel casters
- Genesis casters are maintenance-free (#62-B)



P-1015

# footstools

**features**

- Strong, all welded steel
- Legs and handrail constructed of 1" O.D. 16-gauge tubular steel
- Corrugated, black rubber top surface
- Rubber floor tips



**P-10**

Powder Coat

12"D x 16"W x 8"H

300 lbs. weight capacity



**P-10-A**

Powder Coat

12"D x 16"W x 8"H

300 lbs. weight capacity

- 39"H with handrail



**P-1010-A-SS**

Stainless Steel

12"D x 16"W x 8"H

300 lbs. weight capacity

- 39"H with handrail



**P-1010-SS**

Stainless Steel

12"D x 16"W x 8"H

300 lbs. weight capacity



**P-18**

Powder Coat

12"D x 16"W x 16¼"H

300 lbs. weight capacity

- First step 9"H
- Second step 16¼"H



**P-18-A**

Powder Coat

12"D x 16"W x 16¼"H

300 lbs. weight capacity

- First step 9"H
- Second step 16¼"H
- 47"H with handrail



**P-11**

Powder Coat

12"D x 30"W x 8"H

300 lbs. weight capacity



**P-1011-SS**

Stainless Steel

12"D x 30"W x 5"H

300 lbs. weight capacity



**P-1012-SS**

Stainless Steel

12"D x 30"W x 8"H

300 lbs. weight capacity



**P-1015**

Stainless Steel

14"D x 19"W x 6¼"H

- 650 lbs. weight capacity
- Interlock footstools to create larger platforms
- Heavy gauge legs and handles

Safely stack for increased height