

shown with optional accessories

CDS-270 Cart

Overall Size

 $50^{5}/8$ "W x $26^{1}/8$ "D x 65"H

- Designed to be used as a clean/dirty linen cart or as a large utility cart
- May be ordered with optional shelves, wire backs, and other accessories

multi-purpose carts

features

- All welded stainless steel construction for one single piece frame
- 6" chrome plated casters with diagonally mounted swivel locks provide stability
- Wrap around bumper mounted away from frame for easy cleaning
- Optional shelves and accessories available to configure the cart to suit individual needs







Cart shown with optional accessories

CDS-262 Multi Use Procedure Cart

Overall Size

295/8"W x 295/8"D x 623/8"H

- Designed to be used as a surgical or obstetrical case cart
- Multi-purpose cart can also be used as a linen cart or as a clean/dirty utility cart
- May be ordered with optional shelves, wire backs, and other accessories