

Stainless Steel Twin Basin Stand with Shelf

Model SS8360

UMF Medical's Basin Stands are constructed on heavy duty tubular stainless steel framework. Rubber bumpers inside the basin holders eliminate rattle and scarring of the basins. The units are mounted on swivel casters and come with a 3-year limited warranty.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › Heavy duty stainless steel tubular framework
- › Two 13" (33cm) basins, 7 Quart (6.6L) capacity
- › Rubber bumpers inside basin holder
- › Single shelf below basins
- › 2" (5.1cm) ball bearing swivel casters
- › 3 year limited warranty

Specifications

Height: 14" (35.6cm)
Base Diameter: 14.5" (36.8cm)
Pail capacity: 13 Qt. (12.3L)