



SCC-256-MS shown with optional drop leaf

SCC-256-MS

Overall outside dimensions:
28⁵/₈"W x 28⁵/₈"D x 39⁷/₈"H

Usable Interior Space
20¹/₄"W x 24"D x 26¹/₄"H

One shelf configuration recommended

sealed carts

surgical

features

Sealed Case Carts features a full seal around the door perimeter and center stop providing an added level of protection against dust and debris.

- Fully-welded all stainless steel construction
- Exterior push handles for better control and more interior space
- One-piece continuous wrap-around bumper, non-marring
- Sound-deadened tops and doors for quiet operation
- Positive lock door latches
- All carts can be manufactured to individual specifications



Biohazard Indicator Wheel
CartStat



Sealed Case Cart Option



Heavy Duty Casters
Option



Floor Sloped Toward Rear
2° slope for containment of fluid.



Total Lock Brake
Option



SCC-235 shown with optional shelves & CartStat

SCC-242 shown with optional shelves

SCC-245 shown with optional shelves

SCC-233

Overall outside dimensions:
28⁵/₈"W x 28⁵/₈"D x 56⁷/₈"H

Usable Interior Space
20¹/₄"W x 24"D x 43¹/₄"H

Two shelf configuration recommended

SCC-235

Overall outside dimensions:
34¹/₈"W x 32³/₈"D x 56⁷/₈"H

Usable Interior Space
25³/₄"W x 27³/₄" x 43¹/₄"H

Two shelf configuration recommended

SCC-242

Overall outside dimensions:
45⁷/₈"W x 28⁵/₈"D x 39⁷/₈"H

Usable Interior Space
37¹/₂"W x 24"D x 26¹/₄"H

One shelf configuration recommended

SCC-245

Overall outside dimensions:
45⁷/₈"W x 28⁵/₈"D x 56⁷/₈"H

Usable Interior Space
37¹/₂"W x 24"D x 43¹/₄"H

Two shelf configuration recommended

OPTIONAL FEATURES

- Security Latch
- Heavy Duty 6" x 2" Casters
- Color Identification Wheel
- Barcode Plate Holder
- 8¹/₂" x 11" Card Holder

Shelving Configurations

CDS-233-SS Stainless steel solid shelf for CDS-233, CDS-256-LD and CDS-256-MS.

CDS-233-SSRO ... Stainless steel roll-out solid shelf for CDS-233, CDS-256-LD and CDS-256-MS.

CDS-233-WS Stainless steel wire shelf for CDS-233, CDS-256-LD and CDS-256-MS.

CDS-233-WSRO .. Stainless steel roll-out wire shelf for CDS-233, CDS-256-LD and CDS-256-MS.

CDS-235-SS Stainless steel solid shelf for CDS-235.

CDS-235-SSRO ... Stainless steel roll-out solid shelf for CDS-235.

CDS-235-WS Stainless steel wire shelf for CDS-235.

CDS-235-WSRO .. Stainless steel roll-out wire shelf for CDS-235.

CDS-242-SS Stainless steel solid shelf for CDS-242 and CDS-245.

CDS-242-SSRO ... Stainless steel roll-out solid shelf for CDS-242 and CDS-245.

CDS-242-WS Stainless steel wire shelf for CDS-242 and CDS-245.

CDS-242-WSRO .. Stainless steel roll-out wire shelf for CDS-242 and CDS-245.

Custom SCC carts engineered to customer specifications available upon request.